## ROTOFIX 32 A

Ⓝ DE	Bedienungsanleitung .....	5
Ⓝ EN	Operating Instructions .....	18
Ⓝ FR	Mode d'emploi .....	30
Ⓝ IT	Istruzioni per l'uso .....	43


Fig. 1


Fig. 2 ROTOFIX 32 A

**EG-Konformitätserklärung  
EC Conformity Declaration  
Déclaration de conformité CE  
Dichiarazione di conformità alle norme CEE**

Andreas Hettich GmbH & Co. KG • Föhrenstraße 12 • D-78532 Tuttlingen • Germany

Das bezeichnete Gerät, inklusive Zubehör entspricht den aufgeführten EG-Richtlinien und Normen.  
The denoted device, including accessories corresponds to the listed EC guidelines and standards.  
L'appareil désigné, y compris les accessoires, correspond aux directives CE et aux normes énumérées.  
L'apparecchio designato, compresi gli accessori, è conforme alle direttive CE e alle norme citate.

Geräteart, Type of device, Type d'appareil, Tipo di apparecchio:

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**Contents**

1	Use according to specification .....	19
2	Residual risks .....	19
3	Technical specifications .....	19
4	Notes on safety .....	20
5	Symbol meanings .....	21
6	Delivery checklist .....	21
7	Unpacking the centrifuge .....	21
8	Initial operation .....	21
9	Opening and closing the lid .....	22
9.1	Opening the lid .....	22
9.2	Closing the lid .....	22
10	Installation and removal of the rotor .....	22
11	Loading the rotor .....	22
12	Control and display elements .....	23
12.1	Symbols on the control panel .....	23
12.2	Keys and setting options .....	23
13	Setting the brake step .....	24
14	Setting the centrifuging radius .....	24
15	Centrifugation .....	24
15.1	Centrifugation with preselected time .....	25
15.2	Continuous operation .....	25
15.3	Short-time centrifugation .....	25
15.4	Display of the relative centrifugal force (RCF) .....	25
16	Relative centrifugal force (RCF) .....	25
17	Centrifugation of materials with higher density .....	26
18	Rotor Identification .....	26
19	Emergency release .....	26
20	Maintenance and servicing .....	26
20.1	Centrifuge .....	26
20.2	Rotors and Attachments .....	27
20.2.1	Trunnions .....	27
20.2.2	Rotors and accessories with limited term of use .....	27
20.3	Autoclaving .....	27
20.4	Centrifuge containers .....	27
21	Faults .....	28
22	Change mains input fuse .....	29
23	Acceptance of the centrifuges for repair .....	29
24	Disposal .....	29
25	Anhang / Appendix .....	56
25.1	Rotoren und Zubehör / Rotors and accessories .....	56

### 1 Use according to specification

The machine presented here is a medical product (laboratory centrifuge) according to the IVD guideline 98/79/EG. The centrifuge is used to separate substances or substance mixtures with a density of max. 1.2 kg/dm<sup>3</sup>. This also includes substances and substance mixtures of human origin. The centrifuge is only intended to be used for this purpose. A different use or application over and above this is deemed not in accordance with the specifications. The company Andreas Hettich GmbH & Co. KG undertakes no liability for damages resulting therefrom.

Belonging to the application according to specification is also the observance of all references contained in the Instruction Manual and compliance with the inspection and maintenance works.

### 2 Residual risks

The machine is constructed according to the state of the art and the recognized technical safety regulations. Improper use and handling can result in dangers to life and limb of the user or third parties and impairments to the machine or to other material assets. The machine is only to be used for the specified applications and only in an impeccable technical safety condition.

Disturbances that can interfere with the safety are to be immediately rectified.

### 3 Technical specifications

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen	
Model	ROTOFIX 32 A	
Type	1206	1206-01
Mains voltage (± 10%)	208 – 240 V 1~	100 – 127 V 1~
Mains frequency	50 – 60 Hz	50 – 60 Hz
Connected load	300 VA	300 VA
Current consumption	1.4 A	3.0 A
Max. capacity	4 x 100 ml / 32 x 15 ml	
Allowed density	1.2 kg/dm <sup>3</sup>	
Speed (RPM)	6000	
Force (RCF)	4186	
Kinetic energy	3160 Nm	
Obligatory inspection (BGR 261)	no	
Ambient conditions (EN 61010-1)	<ul style="list-style-type: none"> <li>– Set-up site Indoors only</li> <li>– Altitude Up to 2000 m above sea level</li> <li>– Ambient temperature 2°C to 40°C</li> <li>– Humidity Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C.</li> <li>– Excess-voltage category (IEC 60364-4-443) II</li> <li>– Pollution degree 2</li> </ul>	
Device protection class	I	
Not suitable for use in explosion-endangered areas.		
EMC	<ul style="list-style-type: none"> <li>– Emitted interference (suppression of radio interference) EN 55011, Group 1, Class B EN 61000-3-2 EN 61000-3-3</li> <li>– Interference immunity EN 61000-6-2</li> </ul>	
		FCC Class B
		----
Noise level (dependent on rotor)	≤ 57 dB(A)	
Dimensions	<ul style="list-style-type: none"> <li>– Width 366 mm</li> <li>– Depth 430 mm</li> <li>– Height 257 mm</li> </ul>	
Weight	23 kg	

#### 4 Notes on safety

**No claim under guarantee will be considered by the manufacturer unless the above instructions have been adhered to.**

- **Before the initial operation of your centrifuge you should read and pay attention to the operating instructions.  
Only personnel that has read and understood the operating instructions are allowed to operate the device.**
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner. These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate.
  - However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge should be installed on a good, stable base.
- Before using the centrifuge absolutely check the rotor for firm placement.
- When the centrifuge is running, according to IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either heat up for at least 3 hours in the warm room before being connected to the mains, or run hot for 30 minutes in the cold room.
- Only the rotors and accessories approved by the manufacturer for this device may be used (see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").
- The centrifuge rotor may only be loaded in accordance with the chapter "Loading the rotor".
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm<sup>3</sup>.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with:
  - inflammable or explosive materials
  - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.
 

For hazardous substances centrifuge containers with special screw caps must strictly be used. In addition to the screw cap centrifuge containers, for materials in hazard category 3 and 4 a biosafety system must be used (see the World Health Organisation's "Laboratory Biosafety Manual").

Under a biosafety system small drips and aerosols are prevented from escaping by a bioseal (packing ring) located between the hanger and the lid.

If the hanger of a biosafety system is used without the lid, the packing ring must be removed from the hanger in order to prevent the packing ring from being damaged during the centrifugation run. Damaged packing rings must not be used to seal the biosafety system.

Without the use of a biosafety system the centrifuge is not microbiologically sealed in the sense of the EN 61010-2-020 standard.

For further details of available biosafety systems see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". If in doubt, you should obtain relevant information from the manufacturer.
- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.
- Rotors, suspensions and accessories that possess traces of corrosion or mechanical damage or if their term of use has expired may not be used any longer.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Andreas Hettich GmbH & Co. KG company are allowed to be utilised.
- The following safety regulations apply:
  - IEC 61010-1 and IEC 61010-2-020 as well as their national deviations.
- The safe operation and reliability of the centrifuge can only be guaranteed if:
  - the centrifuge is operated in accordance with the operating instructions,
  - the electrical installation on the site where the centrifuge is installed conforms to the demands of IEC stipulations,
  - prescribed tests to BGV A1, BGR 261 are carried out by an expert.

## 5 Symbol meanings

Symbol on the machine:

Attention, general hazard area.

Before using the centrifuge implicitly read the operating instructions and pay attention to the safety relevant references!

Symbol in the operating instructions:

Attention, general hazard area.

This symbol refers to safety relevant warnings and indicates possibly dangerous situations.

The non-adherence to these warnings can lead to material damage and injury to personal.

Symbol in the operating instructions:

This symbol refers to important circumstances.

Symbol on the machine and in the operating instructions:

Symbol for the separate collection of electric and electronic devices according to the guideline 2002/96/EG (WEEE). The device belongs to Group 8 (medical devices).

Applies in the countries of the European Union, as well as in Norway and Switzerland.

## 6 Delivery checklist

The following items and accessories are delivered with the centrifuge:

- |   |                                      |
|---|--------------------------------------|
| 1 | Connecting cable                     |
| 2 | Fuses                                |
| 1 | Lubricating grease for trunnions     |
| 1 | Hex. pin driver                      |
| 1 | Release pin                          |
| 1 | Notes on moving the equipment safely |
| 1 | Operating instructions               |

The rotor(s) and associated accessories are included in the delivery in the quantity.

## 7 Unpacking the centrifuge

- Lift the carton upward and remove the padding.

Do not lift by the handle rail.

Observe the weight of the centrifuge, refer to chapter "Technical specifications".

Lift the centrifuge on both sides with an appropriate number of helpers and place it on the laboratory table.

## 8 Initial operation

- According to the laboratory instrument standards IEC 61010-2-020 an emergency switch to separate power supply in the event of a failure must be installed in the building electrical system. This switch has to be placed remote from the centrifuge, preferred outside of the room in which the centrifuge is installed or near by the exit of this room.
- Remove the transportation safety device from the bottom of the housing, see sheet "Transportation safety device".
- Position the centrifuge in a stable and level manner in a suitable place. During set-up, the required safety margin of 300 mm around the centrifuge is to be kept according to IEC 61010-2-020.

When the centrifuge is running, according to IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

- Do not place any object in front of the ventiduct. Keep a ventilation area of 300 mm around the ventiduct.
- Check whether the mains voltage tallies with the statement on the type plate.
- Connect the centrifuge with the connection cable to a standard mains socket. For connection ratings refer to Chapter "Technical specifications".
- Turn on the mains switch. Switch position "I". The last used centrifuge data will be displayed.
- Open the lid.

## 9 Opening and closing the lid

### 9.1 Opening the lid

The lid can only be opened when the centrifuge is switched on and the rotor is at rest. If it cannot be opened under these circumstances, see the section on "Emergency release".

- Swing handle rail on the lid upwards. The symbol "L" (lid open) illuminates in the rotation indicator Ⓞ.
- Open the lid.


### 9.2 Closing the lid

Do not bang the lid shut.

- Place the lid and swing handle rail on the lid downward. The symbol "—" (lid closed) illuminates in the rotation indicator Ⓞ.

## 10 Installation and removal of the rotor

- Clean the motor shaft (C) and the rotor drilling (A), and lightly grease the motor shaft afterwards. Dirt particles between the motor shaft and the rotor hinder a perfect seating of the rotor and cause an irregular operation.
- Place the rotor vertically on the motor shaft. The motor shaft dog (D) has to fit in the rotor slot (B). The alignment of the groove is labelled on the rotor.
- Tighten the rotor tension nut with the supplied wrench by turning in a clockwise direction.
- Check the rotor for firm seating.
- Loosening the rotor: Loosen the tension nut by turning in a counter clockwise direction, and turning until the working point for lifting. After passing the working point for lifting the rotor is loosened from the motor shaft cone. Turn the tension nut until the rotor is able to be lifted from the motor shaft.

## 11 Loading the rotor

Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2).

- Check the rotor for firm seating.
- With swing-out rotors all rotor positions must be lined with **identical** hangers. Certain hangers are marked with the number of the rotor position. These hangers may only be used in the respective rotor position.
- The rotors and hangers may only be loaded symmetrically. For authorised combinations see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".  
In the case of angle rotors all possible rotor positions must be loaded, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".
- On certain suspensions, the weight of the maximum load and the maximum weight of the suspension when it is fully equipped is displayed. This weight may not be exceeded. The weight specified for the maximum loading includes the total weight of adapter, frame, centrifuging container and content.
- In containers with rubber inserts, the same number of rubber inserts must always be among the centrifuge containers.
- Always fill the centrifuge containers outside of the centrifuge.
- No liquid should be allowed to enter the centrifugal chamber during filling and swinging out of the hangers.
- The maximum filling quantity for the centrifuge containers specified by the manufacturer must not be exceeded.
- In order to maintain the weight differences within the centrifuge container as marginal as possible, a consistent fill level in the containers is to be heeded.

## 12 Control and display elements

See figure on page 2.

Fig. 2: Display and control panel

### 12.1 Symbols on the control panel

Rotation indicator. The rotation indicator lights up and rotates anticlockwise while the rotor is turning.

When the rotor is stationary, the status of the lid is displayed by symbols in the rotation indicator:

Symbol  : Lid open

Symbol  : Lid closed

Operator errors and occurring faults are indicated on the display (see Chapter "Faults").

### 12.2 Keys and setting options

RPM/RCF x 100

- Speed

A numeric value of 500 RPM up to the maximum rotor speed can be set. For maximum rotor speed, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". Preset in steps of 100 (RPM = displayed value x 100).

If the key  or  is kept pressed, the value changes with increasing speed.

- Display the brake step and the centrifuging radius.

t

- Running time

- Preset from 1 - 99 minutes, in 1 minute steps

- Continuous operation "--"

- Centrifuging radius. Input in centimeters. Preset from 5 - 16 centimeters, in 1 centimeter steps. For centrifuging radius, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".
- Braking steps 0 or 1. Step 1 = short run-down time, Step 0 = long run-down time.

If the key  or  is kept pressed, the value changes with increasing speed.

- Start centrifugation run.

- End centrifugation run.  
The rotor runs down with the preselected brake step.
- Save the brake step and the centrifuging radius.


- Display of the relative centrifugal force (RCF).  
The display of the relative centrifugal force (RCF) appears while the key  is kept pressed.


- Short-time centrifugation.  
The centrifugation run occurs while the key  is kept pressed.
- Display the brake step and the centrifuging radius.

### 13 Setting the brake step

- Switch off the mains switch.
- Keep the key  beneath the speed indicator and the key **[IMPULS]** pressed simultaneously.
- Switch on the mains switch and release the keys again.
- Press the key  beneath the speed indicator until the following display appears:


The machine version set in the factory (e.g. 7) is displayed in the speed indicator and the set brake step is displayed in the time indicator.

- Set the desired brake step with the keys   beneath the time indicator.  
Step 1 = short run-down time, Step 0 = long run-down time.  
For run-down times, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".
- Press the key **[STOP]** to save the setting.


### 14 Setting the centrifuging radius

 The centrifuging radius must be entered in centimeters.

- Switch off the mains switch.
- Keep the key  beneath the speed indicator and the key **[IMPULS]** pressed simultaneously.
- Switch on the mains switch and release the keys again.
- Press the key  beneath the speed indicator until the following display appears:


The set centrifuging radius is displayed in the speed indicator.

- Set the desired centrifuging radius with the keys   beneath the time indicator.  
For centrifuging radius, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".
- Press the key **[STOP]** to save the setting.

### 15 Centrifugation

 When the centrifuge is running, according to IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

 If the permissible weight difference is exceeded within the rotor loading, the drive switches off during the run-up time, and error **-3-** is displayed (see chapter "Faults").  
The centrifugation run can be interrupted at any time by pressing the key **[STOP]**.  
The time and speed can be changed during the centrifugation run, with the keys  .  
If the key  or  is kept pressed, the value changes with increasing speed.  
After a centrifugation run, the display flashes until the cover is opened or a key is pressed.  
If the symbol "—" (lid closed) and "L" (lid open) flashes alternately in the rotation indicator , operation of the centrifuge can only be continued after opening the lid.  
If **rot xx** is displayed, no centrifugation run has taken place because the rotor has been changed, see chapter "Rotor Identification".

- Switch on the mains switch (switch position "I").
- Load the rotor and close the centrifuge cover.

### 15.1 Centrifugation with preselected time

- Set the desired speed with the keys  $\blacktriangle$   $\blacktriangledown$  beneath the speed indicator.
- Set the desired time with the keys  $\blacktriangle$   $\blacktriangledown$  beneath the time indicator.
- Press the key **[START]**. The rotation indicator  $\odot$  appears while the rotor is turning.

The time is displayed in minutes. The last minute is counted down in seconds.  
When the time is displayed in minutes, a point flashes next to the number.

- After expiry of the time or if the centrifugation run is interrupted by pressing the key **[STOP]**, the rotor runs down with the set brake step.

During the centrifugation run, the rotor speed or the resulting RCF value and the remaining time are displayed.

### 15.2 Continuous operation

- Set the desired speed with the keys  $\blacktriangle$   $\blacktriangledown$  beneath the speed indicator.
- Set the time to zero with the key  $\blacktriangledown$  beneath the time indicator. "--" is displayed.
- Press the key **[START]**. The rotation indicator  $\odot$  appears while the rotor is turning. The time count starts from 0.

The first minute is counted up in seconds, and then the time is displayed in minutes.  
When the time is displayed in minutes, a point flashes next to the number.

- Press the key **[STOP]** to end the centrifugation run. The rotor runs down with the set brake step.

During the centrifugation run, the rotor speed or the resulting RCF value and the expired time are displayed.

### 15.3 Short-time centrifugation

- Set the desired speed with the keys  $\blacktriangle$   $\blacktriangledown$  beneath the speed indicator.
- Keep the key **[IMPULS]** pressed. The rotation indicator  $\odot$  appears while the rotor is turning. The time count starts from 0.

The first minute is counted up in seconds, and then the time is displayed in minutes.  
When the time is displayed in minutes, a point flashes next to the number.

- Release the key **[IMPULS]** again to end the centrifugation run. The rotor runs down with the set brake step.

During the centrifugation run, the rotor speed and the expired time are displayed.

### 15.4 Display of the relative centrifugal force (RCF)

The relative centrifugal force (RCF) can be displayed during the centrifugation run.

If the relative centrifugal force (RCF) is used, the centrifuging radius must be entered.

- Keep the key **[RCF]** pressed during the centrifugation run.  
The relative centrifugal force (RCF) appears in the speed indicator (RCF = displayed value x 100).
- Release the key **[RCF]** again. The speed is displayed.

## 16 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left( \frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge.

For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

## 17 Centrifugation of materials with higher density

The rotors are designed to centrifuge substances up to a maximum mean homogenous density of 1.2 kg/dm<sup>3</sup> when rotating at the stated speed.

Denser substances must be centrifuged at lower speed.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density}}} \times \text{Rated speed}$$

e.g.: RPM 4000, density 1.6 kg/dm<sup>3</sup>

$$n_{\text{red}} = \sqrt{\frac{1.2}{1.6}} \times 4000 = 3464 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

## 18 Rotor Identification

After every start of a centrifugation run the rotor utilised is identified.

After a rotor change, the drive switches off and the rotor code (rot xx) is displayed.

- Press the key **START**. The last used centrifuge data will be displayed..

A further operation of the centrifuge is only possible after a single opening of the lid.

If, following a rotor change, the maximum speed of the rotor is less than the set speed, the speed is limited to the maximum speed of the rotor.

## 19 Emergency release

The lid cannot be opened during power failure. An emergency release has to be executed by hand.

For emergency release disconnect the centrifuge from the mains.

Open the lid only during rotor standstill.

Only the plastic release pin provided may be used for emergency release.

See figure on page 2.

- Switch off the mains switch (switch position "0").
- Look through the window in the lid to be sure that the rotor has come to a standstill.
- Insert the release pin (see scope of supply) horizontally into the hole (Fig. 1, A). Push the unlocking pin in until the handle can be lifted when the pin is pressed down.
- Open the lid.

## 20 Maintenance and servicing

Pull the mains plug before cleaning.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Cleaning agents and disinfectants which lie in the pH range 5 – 8 are to be utilised. Alkaline cleaning agents with a pH value > 8 are to be avoided.
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

### 20.1 Centrifuge

- Clean the centrifuge housing and the centrifuging chamber regularly, using soap or a mild detergent and a damp cloth if required. For one thing, this services purposes of hygiene, and it also prevents corrosion through adhering impurities.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately. For surface disinfection we recommend Bacillol® manufactured by Bode Chemie in Hamburg or Biocidal ZF™ from the company WAK-Chemie Medical GmbH in Steinbach.
- Lightly grease the rubber seal of the centrifugal chamber after every cleaning.

## 20.2 Rotors and Attachments

- In order to prevent corrosion and material changes, rotors and accessories must be cleaned regularly with soap or a mild detergent and a damp cloth. Cleaning is recommended at least once a week, even better after every usage.
- If the rotor or accessory parts are contaminated by pathogenic or radioactive material, a suitable cleaning has to be executed. For disinfection we recommend Helipur® H plus N from the company B. Braun Melsungen. For the removal of radioactive material we recommend decon neutracon® from the company Decon Laboratories Limited.
- The rotors and accessory parts must be dried immediately after cleaning.
- Angle rotors, container and hanger made of aluminium are to be lightly greased after drying using acid-free grease, e.g. vaseline.
- With aerosol-sealed rotors and bio safety systems (see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories") the sealing rings are to be checked and cleaned regularly (weekly). The sealing ring is to be replaced immediately upon indication of crack formation, embrittlement or abrasive wear. To prevent the sealing ring from twisting when the lid is opened and closed, the sealing ring must be lightly greased with acid-free grease, such as Vaseline.
- In order to prevent corrosion as a result of moisture between the rotor and the motor shaft, the rotor should be disassembled and cleaned at least once a month, and the motor shaft should be lightly greased.
- The rotors and accessory parts are to be checked on a monthly basis for corrosion damage.

Rotors and attachments may no longer be utilised upon indication of wear and tear or corrosion.

- Check the firm seating of the rotor on a weekly basis.

### 20.2.1 Trunnions

With swing-out rotors the trunnions must be regularly lubricated (Hettich Lubricating Grease No. 4051) in order to ensure consistent swinging out of the hangers.

### 20.2.2 Rotors and accessories with limited term of use

The use of specific rotors, suspensions and accessories is time limited.

They are marked with an expiry date, e.g. "einsetzbar bis Ende: / usable until end of: IV. Quartal 2011" (applicable until the end of: IVth quarter 2011).

The rotors, suspensions and accessories may not be used for longer periods for safety reasons once the marked expiry date has been reached.

## 20.3 Autoclaving

Swing-out rotors, angle rotors made of aluminium, suspension made of metal, lids with biodegradable seals as well as stands and reductions can be autoclaved at 121° C / 250°F (20 mins).

Otherwise you must ask the manufacturer.

The lids of the rotors and containers must be removed prior to autoclaving. Autoclaving accelerates the ageing process of plastics. In addition, autoclaving may discolour plastics. After autoclaving, we recommend that the sealing rings of the aerosol-tight and bio-safety systems be exchanged.

## 20.4 Centrifuge containers

- With leakiness or after the breakage of centrifuging containers broken container parts and leaked centrifugation material are to be completely removed.
- The rubber inserts as well as the plastic sleeves of the rotors are to be replaced after a glass breakage.

Remaining glass splitters cause further glass breakage!

- If this concerns infectious material, a disinfection process is to be executed immediately.

## 21 Faults

If the fault cannot be eliminated with the help of the fault table, please inform Customer Service.

Please state the type of centrifuge and the factory serial number. Both values are visible on the centrifuge type plate.

Perform a MAINS RESET:


- Switch off the mains switch (switch position "0").
- Wait at least 10 seconds and then switch on the mains switch again (switch position "1").

Fault	Display	Cause of fault	Remedy
No display	---	No voltage Mains input fuses defective.	<ul style="list-style-type: none"> <li>- Check distribution voltage.</li> <li>- Check mains power input fuse, refer to Chapter "Change mains input fuse".</li> <li>- Mains switch ON.</li> </ul>
Tacho error	- 1 -	Failure of speed impulses during operation.	<ul style="list-style-type: none"> <li>- Due to safety reasons the lid can only be opened after approx. 120 seconds.</li> <li>- After this time has expired perform a MAINS-RESET.</li> </ul>
System reset	- 2 -	Power failure during the centrifugation run. (The centrifugation run was not finished.)	<ul style="list-style-type: none"> <li>- When stationary, open lid and press <b>[START]</b> key.</li> <li>- Repeat the centrifugation run if necessary.</li> </ul>
Balance error	- 3 -	The rotor is unevenly loaded.	<ul style="list-style-type: none"> <li>- Open lid when rotor is stationary.</li> <li>- Check the loading of the rotor, see chapter "Loading the rotor".</li> <li>- Repeat the centrifugation run.</li> </ul>
Communication	- 4 -	Fault in control unit or power unit.	<ul style="list-style-type: none"> <li>- Perform a MAINS RESET when the rotor has been stationary.</li> </ul>
Overload	- 5 -	Motor or motor control defective.	<ul style="list-style-type: none"> <li>- Perform a MAINS RESET when the rotor has been stationary.</li> </ul>
Overvoltage Undervoltage	- 6 - - 8 -	Supply voltage outside tolerance (see Technical Data).	<ul style="list-style-type: none"> <li>- Perform a MAINS RESET when the rotor has been stationary.</li> <li>- Check supply voltage.</li> </ul>
Overspeed	- 7 -	Fault in the power unit.	<ul style="list-style-type: none"> <li>- Perform a MAINS RESET when the rotor has been stationary</li> </ul>
Excess temperature	- 9 -	Excess temperature switch in motor has triggered.	<ul style="list-style-type: none"> <li>- When rotor is stationary, open lid using emergency unlocking (see Emergency Unlocking chapter).</li> <li>- Allow motor to cool down.</li> </ul>
Version error	A number is displayed in the time indicator.	Incorrect machine version set, control unit jumps into Setting menu.	<ul style="list-style-type: none"> <li>- Set the number 4 using the <b>[▲]</b> <b>[▼]</b> keys underneath the time indicator.</li> <li>- Press the key <b>[STOP]</b> to save the setting.</li> <li>- Perform a MAINS RESET.</li> </ul>
Controller watchdog	- C -	Fault in control unit.	<ul style="list-style-type: none"> <li>- Perform a MAINS RESET when the rotor has been stationary.</li> </ul>
Lid error	- d -	Error in lid locking or lid closure.	<ul style="list-style-type: none"> <li>- Perform a MAINS RESET when the rotor has been stationary.</li> </ul>
Short circuit	- E -	Short circuit in control unit / power unit.	<ul style="list-style-type: none"> <li>- Perform a MAINS RESET when the rotor has been stationary.</li> </ul>
No rotor code	- F -	No rotor recognition at start. No rotor fitted or defective tacho.	<ul style="list-style-type: none"> <li>- Perform a MAINS RESET when the rotor has been stationary.</li> </ul>
New rotor identified	rot...	see section rotor identification.	<ul style="list-style-type: none"> <li>- Press the key <b>[START]</b>.</li> </ul>

## 22 Change mains input fuse

Switch off the mains switch and separate the centrifuge from the mains!

The fuse holder (A) with the mains input fuses is located next to the mains switch.

- Remove the connecting cable from the machine plug socket.
- Press the snap-fit (B) against the fuse holder (A) and remove.
- Exchange defective mains input fuses.

Only use fuses with the rating defined for the type. See the following table.

- Reinsert the fuse holder until the snap-fit clicks shut.
- Reconnect the centrifuge to the mains supply.

Model	Type	Fuse	Order no.
ROTOFIX 32 A	1206	T 3.15 AH/250V	E997
ROTOFIX 32 A	1206-01	T 5 AH/250V	E914

## 23 Acceptance of the centrifuges for repair

If the centrifuge is returned to the manufacturer for repair, it must be decontaminated and cleaned to protect persons, environment and material.

We reserve the right to accept contaminated centrifuges.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

## 24 Disposal

When you are disposing of the device, the respective statutory rules must be observed.

Pursuant to guideline 2002/96/EC (WEEE), all devices supplied after August 13, 2005 may not be disposed as part of domestic waste. The device belongs to group 8 (medical devices) and is categorized in the business-to-business field.

The icon of the crossed-out trash can shows that the device may not be disposed as part of domestic waste.

The waste disposal guidelines of the individual EC countries might vary. If necessary, contact your supplier.

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
1) Sample temp. rise during maximum speed and 1 hour running time

4) can not be centrifugated when plug is attached, Scal. 10µl-300µl, 15ml, 30ml

5) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".




1611		1131					1132		
<b>Ausschwingrotor 8-fach / Swing out rotor 8-times</b>  									
		Röhrchen / tube							
		---	0501	Sarstedt	Sarstedt	Vacutainer	2079	Sarstedt	Sarstedt
									
Kapazität / capacity	ml	5	6	2,6 - 2,9	2,7 - 3	1,6 - 5,0	10	4,0 - 5,5	4 - 7
Maße / dimensions	∅ x L mm	12/13 x 75	12 x 82	13 x 65	11 x 66	13 x 75	17 x 70	15 x 75	16 x 75
Anzahl p. Red./number p. adapter		1	1	1	1	1	1	1	1
Anzahl p. Rotor / number p. rotor		8	8	8	8	8	8	8	8
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000
RZB / RCF		1914	1914	1914	1914	1914	1914	1914	1914
Radius / radius	mm	107	107	107	107	107	107	107	107
 (97%)	sec	22	22	22	22	22	22	22	22
 1	sec	25	25	25	25	25	25	25	25
 0	sec	104	104	104	104	104	104	104	104
Probenerwärmung/Sample temp. rise	K <sup>1)</sup>	6	6	6	6	6	6	6	6



1611		1643				1644			
<b>Ausschwingrotor 8-fach / Swing out rotor 8-times</b>  									
		Röhrchen / tube							
		0578	Vacutainer	---	Sarstedt	0507	Sarstedt	Vacutainer	0518
									
Kapazität / capacity	ml	7	4 - 7	10	4,5 - 5	15	7,5 - 8,2	8,5 - 10	15
Maße / dimensions	∅ x L mm	12 x 100	13 x 100	13 x 100	11 x 92	17 x 100	15 x 92	16 x 100	17 x 100
Anzahl p. Red./number p. adapter		1	1	1	1	1	1	1	1
Anzahl p. Rotor / number p. rotor		8	8	8	8	8	8	8	8
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000
RZB / RCF		2415	2415	2415	2415	2415	2415	2415	2415
Radius / radius	mm	135	135	135	135	135	135	135	135
 (97%)	sec	22	22	22	22	22	22	22	22
 1	sec	25	25	25	25	25	25	25	25
 0	sec	104	104	104	104	104	104	104	104
Probenerwärmung/Sample temp. rise	K <sup>1)</sup>	9	9	9	9	9	9	9	9

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










1) Sample temp. rise during maximum speed and 1 hour running time  
 3) Centrifuge containers of glass will not stand RCF values exceeding 4000



1620A		Reduzierung / adapter						
<b>Winkelrotor 6-fach / Angle rotor 6-times</b>  		SK 63.08						
		Röhrchen / tube						
		0547	0549		0501	Vacutainer	Sarstedt	
								
Kapazität / capacity	ml	85	85	5	6	1,6 - 5	2,6 - 2,9	
Maße / dimensions	Ø x L	mm	38 x 106	38 x 106	12/13 x 75	12 x 82	13 x 75	13 x 65
Anzahl p. Red./number p. adapter		1	1	2	2	2	2	
Anzahl p. Rotor / number p. rotor		6	6	12	12	12	12	
Drehzahl / speed	RPM	6000	6000	6000	6000	6000	6000	
RZB / RCF	<sup>3)</sup>	4186	4186	3783	3783	3783	3783	
Radius / radius	mm	104	104	94	94	94	94	
 (97%)	sec	19	19	19	19	19	19	
 1	sec	22	22	22	22	22	22	
 0	sec	104	104	104	104	104	104	
Probenerwärmung/Sample temp. rise	K <sup>1)</sup>	5	5	5	5	5	5	

1628		1122			1621				
<b>Ausschwingrotor 12-fach / Swing out rotor 12-times</b>   mit / with 1122  mit / with 1621									
		Reduzierung / adapter							
		Röhrchen / tube							
		2079	Sarstedt	Vacutainer	0507	Sarstedt	Vacutainer	0518	
									
Kapazität / capacity	ml	10	4 - 5,5	4 - 7	15	7,5 - 8,2	8,5 - 10	15	
Maße / dimensions	Ø x L	mm	17 x 70	15 x 75	16 x 75	17 x 100	15 x 92	16 x 100	17 x 100
Anzahl p. Gestell/number p. frame		1	1	1	1	1	1	1	
Anzahl p. Rotor / number p. rotor		12	12	12	12	12	12	12	
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	
RZB / RCF		2254	2254	2254	2683	2683	2683	2683	
Radius / radius	mm	126	126	126	150	150	150	150	
 (97%)	sec	22	22	22	22	22	22	22	
 1	sec	25	25	25	25	25	25	25	
 0	sec	110	110	110	110	110	110	110	
Probenerwärmung/Sample temp. rise	K <sup>1)</sup>	12	12	12	16	16	16	16	

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All rights reserved. No part of this publication may be reproduced without the prior written permission of the copyright owner.